



Davecos® Tourniquet

Disposable Tourniquet

Davecos® Tourniquet

Hygienic – Gentle – Simple – Safe



Hygienic

- » One-time usage protects from transmission of germs, avoids cross-contamination and reduces the risk of nosocomial infections such as MRSA
- » Contributes to hygiene standards for personnel
- » No time-consuming cleaning or disinfection
- » One-patient product

Gentle

- » Gentle to skin due to soft fleece
- » No crushing, constriction or pinching
- » Optimum comfort for the patient (3,5 cm wide, 65 cm long)

Simple

- » Type Velcro® micro fastening for optimum fixing and releasing
- » Easily „re-compressable“ if needed

Safe

- » Elastic band provides optimum compression
- » Latex-free
- » Silicon-free
- » DEHP-free

Extra-light weight for emergency services: 2 g / piece

Made in Germany.



REF	Contents	Packaging Unit	PZN
70.521.00.210	Dispenser, 100 pieces	800 pieces	11102732
70.521.00.211	Re-fill bag, 100 pieces	800 pieces	



P. J. Dahlhausen & Co. GmbH
Emil-Hoffmann-Straße 53
50996 Köln · Germany
Phone: +49(2236)39 13-0
Fax: +49(2236)39 13-409
export@dahlhausen.de
www.dahlhausen.de