



Biopsy skin punches **NEW**

Biopsy skin punches NEW



REF 11.900.01.002

The Biopsy skin punch is used in the dermatology for smaller surgeries. It is used to take a skin sample for subsequent dermatopathology evaluation. The skin punch is equipped with a circular stainless steel blade, protection cap and handle.

- **Blade:** stainless steel
- **Handle:** HIPS (high impact polystyrene)
- **Cap:** PE



REF 11.900.01.006

| Diameter skin punch | Size | Colour | PU | REF |
|---------------------|------|--------|-----|---------------|
| 2 mm | 2 | yellow | 600 | 11.900.01.002 |
| 3 mm | 3 | grey | 600 | 11.900.01.003 |
| 4 mm | 4 | blue | 600 | 11.900.01.004 |
| 5 mm | 5 | pink | 600 | 11.900.01.005 |
| 6 mm | 6 | green | 600 | 11.900.01.006 |
| 8 mm | 8 | purple | 600 | 11.900.01.008 |

Seamless Sharp
Blade

Ribbed sure-grip handle
provides comfort and control

Blade size colour
coded handle



DAHLHAUSEN®

P. J. Dahlhausen & Co. GmbH
Emil-Hoffmann-Straße 53
50996 Köln · Germany
Phone: +49 (0) 2236 3913-0
Fax: +49 (0) 2236 3913-409
export@dahlhausen.de
www.dahlhausen.de