

DEHP / BBP / DBP free




## Oxygen Range

Oxygen Cannula  
Oxygen Mask  
Nebulizer  
Oxygen Catheter  
Connecting Tube

# Oxygen Cannula



REF 01.000.01.601



REF 01.000.01.602



REF 01.000.01.750



REF 01.000.01.751

## Oxygen Cannula Soft Tip Straight

210 cm no-crush tubing and standard connector

Description	Sterile	Pack	REF
Neonate	No	50	01.000.01.601
Pediatric	No	50	01.000.01.602
Adult	No	50	01.000.01.750

## Oxygen Cannula Soft Tip Straight

420 cm no-crush tubing and standard connector

Description	Sterile	Pack	REF
Adult	No	50	01.000.00.319

## Oxygen Cannula Soft Tip Curved

210 cm no-crush tubing and standard connector

Description	Sterile	Pack	REF
Adult	No	50	01.000.01.751

All Oxygen Range Products are Latex & DEHP free.

# Oxygen Mask



REF 01.000.01.120



REF 01.000.08.115



REF 01.000.08.110



REF 01.000.08.130

## Aerosol Oxygen Mask

Under-the-chin style, without tubing

Description	Sterile	Pack	REF
Pediatric	No	50	01.000.01.121
Adult	No	50	01.000.01.120

## Oxygen Mask

Under-the-chin style, without tubing

Description	Sterile	Pack	REF
Neonate	No	50	01.000.08.113
Adult	No	50	01.000.08.115

## Oxygen Mask

Under-the-chin style, 210 cm no-crush tubing with 2 standard connectors

Description	Sterile	Pack	REF
Pediatric	No	50	01.000.01.122
Adult	No	50	01.000.08.110

## Oxygen Mask High Concentration

Under-the-chin style, 210 cm no-crush tubing with standard connector, two side valves (without rebreathing) with rebreather bag

Description	Sterile	Pack	REF
Pediatric	No	50	01.000.08.131
Adult	No	50	01.000.08.130

## Oxygen Mask High Concentration

Under-the-chin style, 210 cm no-crush tubing with standard connector, one side valve (partly rebreathing) with rebreather bag

Description	Sterile	Pack	REF
Adult	No	50	01.000.08.132

# Oxygen Mask

## Tracheotomy



REF 01.000.01.118



REF 01.000.01.119

### Full Face Oxygen Mask

6 mm/22 mm connector and elastic band

Description	Sterile	Pack	REF
Adult, 6 mm connector	No	50	01.000.01.116
Adult, 22 mm connector	No	50	01.000.01.117

### Tracheotomy Oxygen Mask

6 mm connector and elastic band

Description	Sterile	Pack	REF
Pediatric/Adult	No	50	01.000.01.118

### Tracheotomy Oxygen Mask

22 mm connector and elastic band

Description	Sterile	Pack	REF
Pediatric/Adult	No	50	01.000.01.119



# Nebulizer



REF 01.000.08.753

## Nebulizer Set (Aerosol)

Under-the-chin style Mask, 210 cm no-crush tubing with standard + universal all fits connector, Nebulizer 6 ml

Description	Sterile	Pack	REF
Pediatric	No	50	01.000.08.752
Adult	No	50	01.000.08.753

## Nebulizer Set (Aerosol)

Under-the-chin style Mask, 210 cm no-crush tubing with 2 standard connectors, Nebulizer 6 ml

Description	Sterile	Pack	REF
Pediatric	No	50	01.000.08.852
Adult	No	50	01.000.08.853



REF 01.000.08.853

## Nebulizer T-piece Kit

T-mouthpiece, 15 cm corrugated tubing, 210 cm no-crush tubing with 2 standard connectors, Nebulizer 6 ml

Description	Sterile	Pack	REF
Adult	No	50	01.000.08.900



REF 01.000.08.900

# Nebulizer



REF 01.000.08.901

## Nebulizer T-piece Kit

T-mouthpiece, 15 cm corrugated tubing, 210 cm no-crush tubing with standard + universal all fits connector, Nebulizer 6 ml

Description	Sterile	Pack	REF
Adult	No	50	01.000.08.901



REF 01.000.08.903

## Nebulizer T-piece Kit

For incorporation in a respiratory circuit; T-piece with both side 22M/15F connector, 210 cm no-crush tubing with standard + universal all fits connector, Nebulizer 6 ml

Description	Sterile	Pack	REF
Adult	No	100	01.000.08.903



REF 01.000.08.904

## Nebulizer therapy system

For incorporation in a respiratory circuit, T-piece with both sides 22F-22M/15F connector with non-return-valve, 210 cm no-crush tubing with standard + universal all fits connector, Nebulizer 6 ml

Description	Sterile	Pack	REF
Adult	No	50	01.000.08.904



# Nebulizer



REF 01.000.09.100



REF 01.000.08.924

## Nebulizer therapy system

For incorporation in a respiratory circuit, T-piece with both sides 22F-22M/15F connectors, with non-return-valve, 210 cm no-crush tubing and 2 standard connectors, Nebulizer 6 ml

Description	Sterile	Pack	REF
Adult	No	50	01.000.08.905

## Nebulizer therapy closed system

With mouthpiece, Y-piece with expiratory valve, filter and inspiratory valve with T-piece with corrugated tube, 210 cm no-crush tubing and 2 standard connectors, Nebulizer 6 ml

Description	Sterile	Pack	REF
Adult (closed system)	No	50	01.000.09.100

## Nebulizer System (Aerosol)

Under-the-chin style Mask, 210 cm no-crush tubing with standard + universal all fits connector, Nebulizer 20 ml

Description	Sterile	Pack	REF
Adult	No	50	01.000.08.924



# Oxygen Catheter



REF 07.061.08.100

## Oxygen Catheter

Distal eye, fixed compress, stepped (male) connector, 40 cm

	Size	Sterile	Pack	REF
	CH 08	Yes	600	07.060.08.100
	CH 10	Yes	600	07.060.10.100
	CH 12	Yes	600	07.060.12.100
	CH 14	Yes	600	07.060.14.100



REF 07.061.08.110

## Oxygen Catheter

Distal eye, fixed compress, funnel (female) connector, 40 cm

	Size	Sterile	Pack	REF
	CH 08	Yes	600	07.061.08.110
	CH 10	Yes	600	07.061.10.110
	CH 12	Yes	600	07.061.12.110
	CH 14	Yes	600	07.061.14.110

## Oxygen Catheter

Distal eye, lateral eyes, atraumatic opening, loose compress, stepped (male) connector, 40 cm

	Size	Sterile	Pack	REF
	CH 08	Yes	600	07.061.08.200
	CH 10	Yes	600	07.061.10.200
	CH 12	Yes	600	07.061.12.200
	CH 14	Yes	600	07.061.14.200

## Oxygen Catheter

Distal eye, 4 lateral eyes, without compress, stepped (male) connector, 40 cm

	Size	Sterile	Pack	REF
	CH 10	Yes	600	07.059.10.250
	CH 12	Yes	600	07.059.12.250

## Oxygen Catheter

Distal eye, 4 lateral eyes, loose compress, funnel connector, 220 cm

	Size	Sterile	Pack	REF
	CH 14	Yes	150	07.061.14.220



# Oxygen Connecting Tube



REF 01.000.02.002

## Oxygen Connecting Tube

No-crush tubing with 2 standard connectors

Length	Sterile	Pack	REF
210 cm	No	50	01.000.02.002
210 cm	Yes	50	01.000.02.003
426 cm	No	50	01.000.02.014
700 cm	No	25	01.000.02.018



REF 01.000.02.014

## Oxygen Connecting Tube

210 cm no-crush tubing with standard + universal all fits connector

Length	Sterile	Pack	REF
210 cm	No	50	01.000.02.099

## Oxygen Connecting Tube

200 cm no-crush tubing, colour green, with funnel/drain connector

Length	Sterile	Pack	REF
200 cm	Yes	50	07.062.14.111



REF 07.062.14.111

## Oxygen Connecting Tube

200 cm no-crush tubing, colour green, with plain/funnel connector

Length	Sterile	Pack	REF
200 cm	Yes	50	07.062.16.110





# Oxygen Connecting Tube



REF 07.062.00.080


## Oxygen Connecting Tube

Colour green, funnel/funnel connector, 180 cm

	Size	Sterile	Pack	REF
	CH 14	Yes	150	07.062.00.080
	CH 16	Yes	150	07.062.16.180

## Oxygen Connecting Tube

Colour green, funnel/funnel connector, 200 cm


	Size	Sterile	Pack	REF
	CH 16	Yes	150	07.062.16.200



REF 40.000.90.300

## Oxygen Connecting Tube

Colour green, plain tubes both sides, 200 cm

	Size	Sterile	Pack	REF
	CH 24	Yes	75	07.062.00.200

## Bubble Tube, transparent (PVC)


Inner diameter of the tube 3,3 mm, 30 m roll (bulk)

	Size	Sterile	Pack	REF
	CH 10	No	1	40.000.90.300

# Connectors



Stepped (Male) Connector



Funnel (Female) Connector



Standard Connector



Universal All Fits Connector



Drain Connector



Blue Latex Connector



P. J. Dahlhausen & Co. GmbH  
Emil-Hoffmann-Straße 53  
50996 Köln · Germany  
Phone: +49 (0) 2236 3913 - 0  
Fax: +49 (0) 2236 3913 - 409  
[export@dahlhausen.de](mailto:export@dahlhausen.de)  
[www.dahlhausen.de](http://www.dahlhausen.de)