



Protective Covers for Ultrasound Probes

Core Products

General Purpose Probe Covers (RADIOLOGY)

DAHLHAUSEN Ultrasound Probe Covers reduces the risk of infectious disease transmission; protects healthcare professionals in every type of ultrasound environment and helps improve clinical outcomes by maximizing the protection available for both patients and medical staff during procedures.

The Dahlhausen comprehensive line of advanced probe cover solutions facilitates the safe and hygienic use of ultrasound transducers. Whatever branch of medicine, wherever ultrasound is used, the Dahlhausen probe covers offers the key to peace of mind when it matters most.

Description	Material	Size	Sterility	Qty/cs	REF
General purpose probe cover	PU	10 x 61 cm	yes	30	42.300.10.062
General purpose probe cover	PU	10 x 122 cm	yes	30	42.300.10.123
General purpose probe cover	PU	10 x 244 cm	yes	30	42.300.10.245
General purpose probe cover	PU	13 x 61 cm	yes	30	42.300.13.062
General purpose probe cover	PU	13 x 122 cm	yes	30	42.300.13.123
General purpose probe cover	PU	13 x 244 cm	yes	30	42.300.13.245
General purpose probe cover	PU	15 x 61 cm	yes	30	42.300.15.062
General purpose probe cover	PU	15 x 122 cm	yes	30	42.201.13.560
General purpose probe cover	PU	15 x 244 cm	yes	30	42.300.15.245

All above mentioned items are single packed in blue wrap together with 2 rubber bands, 2 tape strips and 1 sachet of transmission gel. Further covers are available. Please contact your medical product advisor.



The Dahlhausen General Purpose Probe Covers offer additional sterile protection when performing routine ultrasound and puncturing procedures. The covers are telescopically folded for easy application of transmission gel, as well as easy application of the cover onto the transducer. All covers have a 3D box end.



Transesophageal (TEE) Probe Covers (CARDIOLOGY)



REF 42.202.25.101

The Dahlhausen Trans Esophageal Echocardiography (TEE) covers feature a tapered design to provide a tight fit to most commonly used TEE-transducers. The form-fitting TEE-covers are pre-loaded on the applicator for ease of application onto the TEE-transducer. The TEE-kit contains a bite-guard, providing a safe passage through the patient's mouth, a gel-filled syringe and a twistlock for fixation.

Description	Material	Size	Sterility	Qty/cs	REF
TEE-probe cover kit	PU	25/11/1,000 mm	yes	10	42.202.25.101
TEE-probe cover kit	latex	25/11/1,000 mm	yes	10	42.202.25.201

Further covers are available. Please contact your medical product advisor.

Endocavity Ultrasound Probe Covers (GYN/URO)



REF 42.200.30.201

The Dahlhausen Endocavity Ultrasound Probe Covers provides a germ barrier for transrectal and transvaginal ultrasound procedures and protects both the clinician and the patient.

Description	Material	Size	Sterility	Qty/cs	REF
Endocavity Probe Cover	PU	30 x 200 mm	no	50	42.200.30.201
Endocavity Probe Cover	latex	30 x 300 mm	no	100	42.200.30.303
Endocavity Probe Cover	latex	33 x 195 mm	no	100	42.200.33.197

Further covers are available. Please contact your medical product advisor.

Ultrasound Transducer Covers

With the ever increasing use of Ultrasound in the healthcare setting it is very important to protect your patients and healthcare workers against Hospital Acquired Infections by using the right Ultrasound Transducer Cover.

DAHLHAUSEN offers a full range of disposable protective covers to minimize the risk of CROSS-CONTAMINATION, protect PATIENTS & MEDICAL STAFF while extending the life of your expensive EQUIPMENT.

Our General Purpose probe covers

- Are delivered sterile
- Are telescopically folded for ease application
- Have a 3D box-end
- Available in a wide range of sizes

Our TEE probe covers

- Are delivered sterile
- Have a tapered design to provide a tight fit
- Are pre-loaded on an applicator, therefore ready to use
- Come with a bite-guard to protect the probe

Our Endocavity probe covers

- Are delivered non-sterile
- Fit most commonly used transvaginal/transrectal probes
- Are packed individually

Dahlhausen Ultrasound Probe Covers provides

- Protection against HAI for both patients and staff
- A safe and hygienic use of ultrasound probes
- A longer life for ultrasound probes



P. J. Dahlhausen & Co. GmbH
Emil-Hoffmann-Straße 53
50996 Köln · Germany
Phone: +49 (0) 2236 3913-0
Fax: +49 (0) 2236 3913-409
export@dahlhausen.de
www.dahlhausen.de