

**NEW**



## Secco™ | Secco Protect™

Faecal management systems

Faecal collection bag

Irrigation sets

Contains faecal contamination and promotes a healthy patient environment

**Secco™ FMS is a closed bowel catheter system for the hygienic drainage of semi-liquid and liquid bowel excretions**

- Reduces skin irritation and pressure ulcers
- Protects patients and staff from infections
- Reduces complications and shortens inpatient stay
- Protects wounds from contamination
- Controls odour and enhances patient comfort
- Reduces the costs caused by faecal incontinence
- Protects patient's dignity




Attachment plate with fastening strap




Disposable syringe, 45 ml



Irrigation (blue) and retention connectors (white)



Sampling port



Catheter localisation marking

## Easy and safe to apply

Secco™ features unique superabsorbent technology to bind infectious liquid faeces. It prevents the spread of infectious material when collecting and disposing of stool. The unique EasyFlush surface prevents clogging of the tube and effectively retains odours.

### Easy and safe to apply

- Uncomplicated handling
- Simple introduction
- Safe disposal



Odour-inhibiting filters on the back of the bag



Low pressure retention cuff



1,500 ml collection bag with superabsorbent pad

## Secco Protect™

The logical further development of the proven stool drainage system

**Secco Protect™ fulfils the well-known features of our long-proven stool drainage system:**

- Reduces skin irritation and pressure ulcers
- Protects patients and staff from infections
- Reduces complications and shortens inpatient stay
- Protects wounds from contamination
- Controls odour and enhances patient comfort
- Reduces the costs caused by faecal incontinence
- Protects patient's dignity



Low pressure retention cuff and catheter localisation marking



Finger insert/insertion aid



Blue: irrigation connector  
White: connector for retention cuff with integrated pressure relief valve



Sampling port



Disposable syringe, 45 ml  
Interchangeable seal for attachment plate

## Secco Protect™

The next generation with even more safety and improved application

### Easy and safe to apply

- Even easier insertion due to the new finger insert on the retention cuff
- The new pressure relief valve prevents overfilling of the retention cuff and offers the patient additional assurance
- More hygienic change of the collection bag due to enclosed interchangeable seals

Secco Protect™ features unique superabsorbent technology to bind infectious liquid faeces. It prevents the spread of infectious material when collecting and disposing of stool. The unique EasyFlush surface prevents clogging of the tube and effectively retains odours.



Odour-inhibiting filters on the back of the bag



1,500 ml collection bag with superabsorbent pad



Attachment plate with interchangeable seal and fastening strap



# Secco™ | Secco Protect™

Delivers more than merely preserving patient dignity



REF 52.000.00.150



REF 52.000.00.100

Faecal incontinence is an unpleasant problem for affected patients and has great importance in maintaining patient dignity, preventing cross-infections and preventing skin lesions.

### Secco™ | Secco Protect™ deliver more than merely preserving patient dignity

- A prolonged hospital stay due to C. difficile results in high costs
- Patients suffering from faecal incontinence have a threefold increased risk of developing a pressure ulcer
- With Secco™ you have a high-quality product group at your disposal

Description	REF
Secco™ – closed stool drainage system	52.000.00.100
Secco Protect™ – closed stool drainage system	52.000.00.150
Replacement bags for closed stool drainage systems, 10 pieces	52.000.00.200

### Set contents

- Flexible bowel catheter with low-pressure cuff (silicone) and irrigation connector (length 170 cm)
- 3 x 1,500 ml collection bags with superabsorbent material
- 3 x Interchangeable seals (for Secco Protect™)
- Attachment plate with velcro straps for fastening to the patient's bed
- Instructions for Use
- 45 ml syringe
- Clamp

# Faecal collection bag



Faecal collection bag, self-adhesive, with connector for drainage bag



Connector for drainage bag

Description	REF
Faecal collection bag, self-adhesive, with connector for drainage bag	52.000.00.201

## Your benefits

- Simple and hygienic solution for patients with diarrhoea
- Secure attachment due to trimmable, self-adhesive skin protection plate, which seals the anal area and protects the skin
- Available with two seals
  - » For thin stools with connector for a secretion bag
  - » For soft stool with resealable opening for cleaning

# Irrigation sets



Irrigation set closed, with valve, sterile

Description	REF
Irrigation set with valve, sterile	05.000.22.665
Irrigation flushing system, sterile	05.000.22.899
Irrigation bag 2,000 ml, sterile	05.000.22.898



Irrigation flushing system, sterile



Irrigation bag 2,000 ml, sterile



P. J. Dahlhausen & Co. GmbH  
Emil-Hoffmann-Straße 53  
50996 Köln · Germany  
Phone: +49(2236)39 13-0  
Fax: +49(2236)39 13-109  
[export@dahlhausen.de](mailto:export@dahlhausen.de)  
[www.dahlhausen.de](http://www.dahlhausen.de)