

Best Allround Treatment for Hospitals



ATRAUMA

Atraumatic suction catheters

 DAHLHAUSEN®

ATRAUMA | ATRAUMATIC SUCTION CATHETERS

The atraumatic suction catheter with cuff and opening with 2 lateral eyes minimises mucous membrane lesions during the suction process. Suction catheters consist of a tube with a defined tip geometry and a connector with funnel.



REF. 07.010.16.143

STRAIGHT CATHETER TIP

	SIZE	LENGTH	PU	PRN	REF
	CH 10	60 cm	100	3656155	07.010.10.143
	CH 12	60 cm	100	3656244	07.010.12.143
	CH 14	60 cm	100	3656273	07.010.14.143
	CH 16	60 cm	100	3656310	07.010.16.143
	CH 10	60 cm	800	10518258	07.011.10.143
	CH 12	60 cm	800	10518264	07.011.12.143
	CH 14	60 cm	800	10518270	07.011.14.143
	CH 16	60 cm	800	10518287	07.011.16.143


HiMi: 01.99.01.007



Endotracheal suction tube

OPTIONAL ACCESSORIES

DESCRIPTION	PU	REF
Endotracheal suction tube - 200 cm with fingertip and 2 clamps, smooth end	30 pieces	05.000.22.657
Endotracheal suction tube - 200 cm with fingertip and 2 clamps, with tube piece	30 pieces	05.000.22.658



Endotracheal suction set for collecting bronchial secretions

DESCRIPTION	PU	PRN	REF
Endotracheal suction set with funnel and fingertip	60 pieces	97487212	07.098.00.000
Endotracheal suction set with funnel and graduated connector	60 pieces	08825042	07.098.00.003
Secretion collector (10 ml)	60 pieces	10033851	07.098.00.001