



Regional anaesthesia

Spinal anaesthesia

Epidural anaesthesia

Quincke cannulas

Spinal anaesthesia



Spinal anaesthesia

When the surgeon August Bier and his assistant Hildebrandt undertook their first attempts at spinal anaesthesia in 1898 they could hardly have imagined the development this form of anaesthesia would take. Today spinal anaesthesia has become established as a largely safe and reliable form of regional anaesthesia. Advancements in local anaesthetics and the cannulas used have also led to this success. DAHLHAUSEN offers users an extensive range of cannulas and sets.



Spinal needles with Quincke cannulas

	ø	Length	Introducer	Pack	REF
	18 G	3,00 inch/76 mm	No	25 pcs.	76.111.18.030
	18 G	3,50 inch/89 mm	No	25 pcs.	76.111.18.035
	19 G	3,00 inch/76 mm	No	25 pcs.	76.111.19.030
	19 G	3,50 inch/89 mm	No	25 pcs.	76.111.19.035
	20 G	1,50 inch/38 mm	No	25 pcs.	76.111.20.016
	20 G	3,00 inch/76 mm	No	25 pcs.	76.111.20.030
	20 G	3,50 inch/89 mm	No	25 pcs.	76.111.20.035
	21 G	3,00 inch/76 mm	No	25 pcs.	76.111.21.030
	21 G	3,50 inch/89 mm	No	25 pcs.	76.111.21.035
	22 G	1,50 inch/38 mm	No	25 pcs.	76.111.22.015
	22 G	3,00 inch/76 mm	No	25 pcs.	76.111.22.030
	22 G	3,50 inch/89 mm	No	25 pcs.	76.111.22.035
	22 G	4,72 inch/120 mm	No	25 pcs.	76.111.22.120
	23 G	3,50 inch/89 mm	No	25 pcs.	76.111.23.035
	25 G	3,00 inch/76 mm	No	25 pcs.	76.111.25.030
	25 G	3,50 inch/89 mm	No	25 pcs.	76.111.25.035
	25 G	3,50 inch/89 mm	Yes	25 pcs.	76.111.25.036
	25 G	4,72 inch/120 mm	Yes	25 pcs.	76.111.25.120
	26 G	3,50 inch/89 mm	Yes	25 pcs.	76.111.26.020
	27 G	3,50 inch/89 mm	Yes	25 pcs.	76.111.27.035
	27 G	4,72 inch/120 mm	Yes	25 pcs.	76.111.27.120

We would be pleased to compile a spinal set according to your requirements.

Pencilpoint

Pencil point spinal needles

Pencil point spinal needles

Our pencil point cannulas offers you all the advantages of an atraumatic spinal cannula. The use of atraumatic cannulas helps largely avoid:

- » postspinal headaches
- » postspinal backaches
- » neurological damage
- » nausea and vomiting

Thanks to a special manufacturing process, our pencil point cannulas offer you faster return flow properties to promptly indicate successful cannulation. We offer you our atraumatic pencil point spinal cannulas in the following designs:



Pencil point cannulas

ø	Length	Introducer	Pack	REF
22 G	3,50 inch/89 mm	No	25 pcs.	76.222.22.035
22 G	3,50 inch/89 mm	Yes	25 pcs.	76.222.22.036
22 G	4,72 inch/120 mm	No	25 pcs.	76.222.22.120
24 G	3,50 inch/89 mm	Yes	25 pcs.	76.222.24.035
24 G	4,72 inch/120 mm	Yes	25 pcs.	76.222.24.120
25 G	3,50 inch/89 mm	Yes	25 pcs.	76.222.25.035
25 G	4,72 inch/120 mm	Yes	25 pcs.	76.222.25.120
26 G	3,50 inch/89 mm	Yes	25 pcs.	76.222.26.035
26 G	4,72 inch/120 mm	Yes	25 pcs.	76.222.26.120
27 G	3,50 inch/89 mm	Yes	25 pcs.	76.222.27.035
27 G	4,72 inch/120 mm	Yes	25 pcs.	76.222.27.120
29 G	3,50 inch/89 mm	Yes	25 pcs.	76.222.29.035

We would also be pleased to offer you individually compiled sets here, too.

Epidural anaesthesia and Sets



Epidural catheter set

consisting of:

catheter, epidural cannula, filter, adapter and LOR syringe

Ø	Pack	REF
16 G	10 pcs.	76.260.44.600
17 G	10 pcs.	76.260.44.700
18 G	10 pcs.	76.260.44.800



Spinal anaesthesia set

consisting of:

Tray, graduated, 2-part	1 wound dressing 7 x 5 cm
5 swabs, plum-size	1 spinal cannula
1 syringe luer, 2 ml, 2-part	25 G x 3 1/2 inch
1 syringe luer, 5 ml, 2-part	with guide cannula
1 cannula no. 1	1 perforated towel 50 x 60 cm
1 cannula no. 12	1 pair sterile gloves
1 clamp	

Set	Pack	REF
Type 153	50 pcs.	94.000.00.253

Do you need a set? Contact us.

We would be pleased to offer you an individual set solution.



P. J. Dahlhausen & Co. GmbH
Emil-Hoffmann-Straße 53
50996 Köln · Germany
Phone: +49(2236)39 13-0
Fax: +49(2236)39 13-109
export@dahlhausen.de
www.dahlhausen.de