

Sterile light handle covers **NEW**

Co-developed and recommended by

Dr. Mach
medical lighting + technology



REF 42.700.01.001



REF 42.700.01.001



REF 42.700.01.002



REF 42.700.01.002

The covers of 12 cm and 15 cm are used to treat the surgical handle in a sterile manner to ensure optimum adjustment of the light source for operation.

Two options are available:

- Single-stage toothed cover (42.700.01.001)
- 3 floors of different toothing-length to enable a slip-proof operation of the lamp handle (42.700.01.002)

Benefits:

- Time and cost saving alternative to the conventional treatment of lamp handles.
- Latex-free disposable product
- Easy handling
- Fits all common surgical light handles

Material: PE

Sterile: Yes

Length: 12 cm/15 cm

PU: 200 pcs.

Length	Box Qty/PU	REF
12 cm	50/200	42.700.01.001
15 cm	50/200	42.700.01.002