



Protective Covers for Ultrasound Probes

Radiology/Anesthesia

Radiology/Anesthesia

DAHLHAUSEN offers a comprehensive portfolio of disposable protective covers and accessories fitting most commonly used ultrasound transducers. Sonic heads are protected by covers. The Dahlhausen ultrasound probe covers provides protection against cross-contamination for both patient and staff, increases efficiency and extends the shelf life of your expensive equipment.



REF 42.300.10.062

The Dahlhausen General Purpose ultrasound probe covers are designed to fit linear, sector and convex transducers. All featuring a soft poly-urethane material and a 3D-tip ensuring no distortion of the ultrasound image. Available sterile and non-sterile. Each protective cover is accompanied by a locking clip which prevents ultrasound gel from leaking out.



REF 42.400.14.061

The new insert-card enables a simple application of the probe cover onto the transducer.



REF 42.400.14.123



Protective covers for sonic heads

Description	Material	Sterility	PU	REF
Protective cover for ultrasound probe, 100 x 610 mm	PU	yes	240	42.300.10.061
Protective cover for ultrasound probe, 100 x 610 mm, with sterile gel	PU	yes	240	42.300.10.062
Protective cover for ultrasound probe, 100 x 1.220 mm	PU	yes	240	42.300.10.122
Protective cover for ultrasound probe, 100 x 1.220 mm, with sterile gel	PU	yes	240	42.300.10.123
Protective cover for ultrasound probe, 100 x 2.440 mm	PU	yes	240	42.300.10.244
Protective cover for ultrasound probe, 100 x 2.440 mm, with sterile gel	PU	yes	240	42.300.10.245
Protective cover for ultrasound probe, 130 x 300 mm	PU	yes	240	42.300.13.030
Protective cover for ultrasound probe, 130 x 300 mm, with sterile gel	PU	yes	240	42.300.13.031
Protective cover for ultrasound probe, 130 x 610 mm	PU	yes	240	42.300.13.061
Protective cover for ultrasound probe, 130 x 610 mm, with sterile gel	PU	yes	240	42.300.13.062
Protective cover for ultrasound probe, 130 x 1.220 mm	PU	yes	240	42.300.13.122
Protective cover for ultrasound probe, 130 x 1.220 mm, with sterile gel	PU	yes	240	42.300.13.123
Protective cover for ultrasound probe, 130 x 2.440 mm	PU	yes	240	42.300.13.244
Protective cover for ultrasound probe, 130 x 2.440 mm, with sterile gel	PU	yes	240	42.300.13.245
Protective cover for ultrasound probe, 150 x 610 mm, with sterile gel	PU	yes	240	42.300.15.062
Protective cover for ultrasound probe, 150 x 1.220 mm, with sterile gel	PU	yes	240	42.201.13.560
Protective cover for ultrasound probe, 150 x 2.440 mm, with sterile gel	PU	yes	240	42.300.15.245
Protective cover for ultrasound probe, 140 x 610 mm, with application aid and sterile gel	PU	yes	240	42.400.14.061
Protective cover for ultrasound probe, 140 x 1.220 mm, with application aid and sterile gel	PU	yes	240	42.400.14.123

Further covers are available. Please contact your medical product advisor.

Accessories

Description	PU	REF
Syringe with ultrasound gel, 10 ml, sterile	200	42.200.00.001
Syringe with ultrasound gel, 10 ml, non-sterile	500	42.200.00.101
Ultrasound gel, 20 ml, sterile	800	42.200.00.010
Ultrasound gel, transparent, 250 ml bottle	24	93.050.00.425
Ultrasound gel, transparent, 500 ml bottle	24	93.050.00.435
Ultrasound gel, transparent, viscous, 500 ml bottle	24	93.050.00.436
Ultrasound gel, transparent, 1 L bottle	12	93.050.00.443
Ultrasound gel, transparent, viscous, 1 L bottle	12	93.050.00.444
Ultrasound gel, transparent, 5 L canister	4	93.050.00.048
Ultrasound gel, viscous, 5 L canister	4	93.050.00.049
Pump for ultrasound gel, 5 L canister	10	93.050.00.449
Ultrasound gel, 5 L cubitainer	1	93.050.00.050
Discharge tap for 5 L cubitainer	30	93.050.00.051
Ultrasound gel, 5 L cubitainer + discharge tap	1	93.050.00.052

Temperature Protection Covers

Description	PU	REF
Protective cover for thermal sensor/temperature probe, latex, 5 x 200 mm, sterile	300	42.200.05.201
Protective cover for thermal sensor/temperature probe, latex, 5 x 200 mm, non-sterile	4.000	42.200.05.200
Protective cover for thermal sensor/temperature probe, latex, 8 x 200 mm, sterile	1.000	42.200.08.401
Protective cover for thermal sensor/temperature probe, latex, 8 x 200 mm, non-sterile	1.000	42.200.08.400
Protective bag to cover devices, 750 x 275 x 1.500 mm, non-sterile, 50 pieces per roll	1	42.700.00.150



P. J. Dahlhausen & Co. GmbH
 Emil-Hoffmann-Straße 53
 50996 Köln · Germany
 Phone: +49 (0) 2236 3913-0
 Fax: +49 (0) 2236 3913-409
export@dahlhausen.de
www.dahlhausen.de