



- Sanadep®**
- » Lint free alcohol swabs, size 4 x 5 cm, soaked in 70 % Ethanol with 0.5 % Chlorhexidine Digluconate
 - » Packed in a handy dispenser with an adhesive pad on the bottom of the container securing a solid placement on the table
 - » Proven efficacy for cleaning and disinfecting the skin prior to injection or blood sampling
 - » Packed by either 150 or 250 swabs per container

Instructions for use

Alcohol swabs for cleansing and disinfecting the intact skin or any other surface.
Make sure the skin or surface remains moist for the minimum exposure time of 30 seconds.
Let air dry.

REF	Description	VSE
78.000.13.100	Sanadep® swabs for skin cleansing	150 pcs.
78.000.13.101	Sanadep® swabs for skin cleansing	250 pcs.