

NEW



Bronchoscopy without disconnection now also possible with OptiFlo 48h and 72h!

Benefits of the OptiFlo adapter (REF. 43.004.60.097) at a glance

- bronchoscopy with the OptiFlo 48h and 72h can now also be performed on the ventilated patient using the new OptiFlo adapter
- no tube disconnection necessary
- continuous ventilation
- simple to use with easy insertion
- no loss of PEEP
- improved PEEP seal
- no saturation loss



Adapter for the OptiFlo 72h




Simple and safe use (above, OptiFlo 48h)



Bronchoscopy (above, OptiFlo 48h)

OptiFlo System

The new generation of closed suction systems



Endotracheal suction tube during use

Versatile system

The OptiFlo closed suction system is equipped for virtually all indications. The OptiFlo comes in the widest range of sizes and lengths with different insertion periods and special features.

Closed-circuit suctioning

Various problems occur during long-term ventilation in the day-to-day clinical setting. We have designed the OptiFlo closed-circuit suctioning system to ensure ventilation is particularly safe and hygienic for patients and caregivers.



OptiFlo 48h

OptiFlo 48h

- Reduced costs thanks to 48-hour insertion period
- Double swivel connector
- Atraumatic catheter tip
- Secretion window
- Extra-soft sleeve
- Bronchoscopy possible without disconnection (with adapter)
- Dosing aerosol adapter



OptiFlo 72h safeclean

OptiFlo 72h safeclean

- Reduced costs thanks to 72-hour insertion period
- Flushing solution does not feed back into the trachea
- Greater hygiene & safety
- Sealable flushing chamber
- Atraumatic catheter tip
- Secretion inspection glass
- Extra-soft sleeve
- Bronchoscopy possible without disconnection (with adapter)
- Dosing aerosol adapter



OptiFlo 168h multilink

OptiFlo 168h multilink

- Reduced costs thanks to 168-hour insertion period
- Flushing solution does not feed back into the trachea
- Greater hygiene & safety
- Sealable flushing chamber
- Atraumatic catheter tip
- Extra-soft sleeve
- Dosing aerosol adapter
- Bronchoscopy access included